

# First stranded cDNA Synthesis Kit (SAMscript)

Cat. No.: GRS01-M

Size: 20 reactions

Store at: -20°C

**For research use only**

## **Description:**

The GeneMark **First stranded cDNA Synthesis Kit (SAMscript)** is provided with all components required to perform first-strand cDNA synthesis. SAMscript Reverse Transcriptase (SAMscript RTase) is genetically engineered version of Moloney Murine Leukemia Virus Reverse Transcriptase (M-MLV RTase) which is a frequent choice for cDNA synthesis because of its ease of use and low intrinsic RNase H activity. In the cDNA synthesis step, RNA was reverse transcribed by SAMscript RTase to produce its cDNA up to 5 kb. The RNase Inhibitor supplied with the kit can protect RNA from degradation.

## **Components:**

This kit is suitable to perform **20** reactions of first-strand cDNA synthesis.

Contents	Volume	Cap
1. SAMscript RTase (4000 U), 200U/μl	20 μl	<b>purple</b>
2. 10X SAMscript RT Buffer	200 μl	<b>purple</b>
3. 2 mM dNTP mixture	200 μl	<b>white</b>
4. 25 μM Oligo (dT)18 Primer	50 μl	<b>green</b>
5. 50 μM Random Primer	100 μl	<b>yellow</b>
6. RNase Inhibitor (40 units/μl)	10 μl	<b>clear</b>
7. Nuclease-free H <sub>2</sub> O	1.8 ml	<b>clear</b>

## General Procedure:

- 1) **Sample denaturation:** Prepare the reaction mix by combining the indicated volumes as followed, then heat RNA/primer mixture at 70°C for 5 min, and then cool immediately on ice for at least 3 min and store on ice.

Nuclease-free H <sub>2</sub> O	X µl
RNA template (0.01 pg-5 µg)	Y µl
Primer *	1 µl
<hr/>	
Total	14.5 µl

\* Template RNA: 10 ng-5 µg Total RNA, 1 ng-0.5 µg mRNA, 0.01 pg-0.5 µg Specific RNA

\*\*Suggestion: each 10 µM gene specific primer; 25 µM Oligo dT Primer; 50 µM Random Primer

- 2) **cDNA synthesis:** Prepare the following reaction mix.

RNA/Primer mix on ice	14.5 µl
10X SAMscript RT Buffer	2 µl
SAMscript RTase	1 µl
2 mM dNTP Mix	2 µl
RNase Inhibitor (40 units/µl)	0.5 µl
<hr/>	
Total	20 µl

- 3) **Incubation condition of the cDNA synthesis reaction:** Based on the primer used in the reaction, use the incubation temperatures and times as below.

**Oligo dT:** 37~42°C for 30~60 min.

**Gene-Specific primer:** 37~42°C for 30~60 min.

**Random Hexamers:** 25°C for 10 min, followed by 37~42°C for 30~60 min.

- 4) **Terminate the Reaction:** Terminate the reaction by incubating at 72°C for 10 min.

\* The cDNA can be directly used in subsequent applications, or stored at -20°C for future use.



**GMbiolab Co., Ltd.**

[www.genemarkbio.com](http://www.genemarkbio.com)

E-mail: [tech@genemarkbio.com](mailto:tech@genemarkbio.com)