



## PRODUCT INFORMATION

**Product Name:** BSA Standart  
**Catalog No.:** GD0069-2, GD0069-10  
**Size:** 2 ml, 10 ml  
**Concentration:** 1mg/ml  
**Mol wt:** ~66 kDa  
**Store at:** -20°C

### **Description:**

BSA Standards are high-quality and universally accepted reference samples for generating accurate standard curves and calibration controls in total protein assays. Useful as an electrophoresis marker and standard for protein concentration reference standards for use in BCA, Bradford and other protein assay protocols.

