



PRODUCT INFORMATION

Product Name:	10X TG Buffer
Catalog No.:	GB04-1 (1 L), GB04-4 (4 L)
Size:	1 L, 4 L
Storage:	RT
Composition:	0.25 M Tris, 1.92 M Glycine, pH 8.3±0.05
Description:	TG buffer is the electrophoresis buffer for polyacrylamide gels, and does not include SDS for protein denature in PAGE. TG Buffer is also use for transfer of Western blot. The buffer is prepared as a 10X concentrate and stored at room temperature.