



PRODUCT INFORMATION

| | |
|----------------------|--|
| Product Name: | 10X PBS Buffer |
| Catalog No.: | GB07-1 (1 L), GB07-4 (4 L) |
| Size: | 1 L, 4 L |
| Storage: | RT |
| Composition: | 1.37 M Sodium Chloride, 27 mM Potassium Chloride, 43 mM $\text{Na}_2\text{HPO}_4 \cdot 7\text{H}_2\text{O}$, 14 mM KH_2PO_4 , pH 7.4 ± 0.05 |
| Description: | PBS buffer is commonly used in biological researches such as cell culturing. The buffer is prepared as a 10X concentrate and stored at room temperature. |