



# PRODUCT INFORMATION

For research use only

---

<b>Product Name:</b>	<b>10X TBST Buffer</b>
<b>Catalog No.:</b>	GB08T-1 (1 L)
<b>Size:</b>	1 L
<b>Storage:</b>	RT
<b>Composition:</b>	0.25M Tris, 1.37M Sodium Chloride, 0.027M Potassium Chloride, 1% Tween 20
<b>Description:</b>	TBST buffer is prepared as a 10X concentrate and stored at room temperature, and is used in Western Blot for protein detection.